

EarPopper™



PROFESSIONAL MODEL

Designed for treatment in the office



EP-3000
EarPopper™
Office Version
Includes 10 disposable tips
One 9-volt battery

EP-1000
Disposable Tips
Non-sterile
Single Use
50 Per Box

HOME VERSION

Designed for continued treatment at home



EP-2100
EarPopper™
Home Version
Includes storage case
One reusable tip
(4) AAA batteries

The EarPopper™ is a safe, simple, non-surgical, non-drug related device for treating such common conditions as Otitis Media with Effusion, Aerotitis/Barotitis (caused by rapid elevation changes), and Eustachian Tube Dysfunction in children and adults.

Based on the Politzer Maneuver, the EarPopper™ delivers a safe, constant, regulated stream of air into the nasal cavity. During the moment of swallowing the air is diverted up the Eustachian tube clearing and ventilating the middle ear. The EarPopper™ relieves negative ear pressure and allows any accumulated fluids to drain.

Designed for treatment in the office with the professional model **EP-3000** and continued treatment at home **EP-2100** Home Version, the EarPopper™ is the safe, simple solution with superior results.



Tecfen Corporation

5860-C Hollister Avenue
Santa Barbara, CA 93117 U.S.A



Tel: (805) 967 1153
Fax: (805) 967 1295



Home: www.tecfe.com
Email: tecfe@tecfe.com