

Ear Wax Removal Kit



60204 Ear Wax Removal Kit

For occasional use as an aid to soften, loosen, and remove excessive ear wax.

Contents of kit:

Ear drops, 15ml bottle
1 pair "EAR PUTTY" Silicone Ear Plugs
Ear Syringe

Safe, Fast-Acting Formula



Packaged 48 kits per Case



Tecfen Corporation

5860-C Hollister Avenue
Santa Barbara, CA 93117 U.S.A



Tel: (805) 967 1153
Fax: (805) 967 1295


Directions for use:



TILT HEAD
SQUEEZE BOTTLE



ROLL "EAR PUTTY" INTO
A BALL-PRESS LIGHTLY



GENTLY FLUSH THE EAR WITH
WARM WATER - USING
EAR SYRINGE



Made in USA

Home: www.tecfe.com

Email: tecfe@tecfe.com