



## **GENERAL SURGICAL**

### TABLE OF CONTENTS

- ◆ Infection Control Glasses, Goggles, Face Masks, and Face Shields
- ◆ General & Utility Surgery Drapes
- ◆ Surgical Skin Marker



# Infection Control Glasses, Goggles, Face Masks, & Face Shields



- Worn Alone
- Vented Sides
- Latex Safe
- Fits over prescription glasses
- Anti-Fog
- 99.9% UV protection



One size

Clear, lightweight plastic glasses provide wrap-around fluid protection and an extremely comfortable fit. Embrace Protective Glasses

Model No. 11000-300

Frame Clear



- UV400 Polycarbonate Lens
- Scratch Resistant Coating
- Lightweight Frames

Model No. 60S43CHS

Frame Clear



Complete peripheral protection with enhanced comfort.

SILO Wear Alone

One size



- Worn Alone
- Vented Sides
- Latex Safe
- Easily adjustable elastic strap
- Optical quality lenses for clarity of vision
- 99.9% UV protection

Sturdy, comfortable, snug-fitting for excellent fluid protection. Vented for additional comfort and reduced fogging.

Model No. 11000-328

Frame Clear

One size


Embrace Protective Goggles



5860-C Hollister Avenue  
Santa Barbara, CA 93117 U.S.A

Tel: (805) 967 1153  
Fax: (805) 967 1295

Home: [www.tecfe.com](http://www.tecfe.com)  
Email: [tecfe@tecfe.com](mailto:tecfe@tecfe.com)



- Worn Alone
- Wrap-around frame
- Secure & comfortable fit
- Non perforated lenses
- Anti-fog on both sides
- Top-down fluid protection ridge
- Fits over prescription glasses
- Latex Safe

The shield can be replaced in a matter of seconds providing extremely economical eye fluid protection.

**Model No.**  
**11000-500**

**Frame**  
**Assorted colors**

**One size**

**Embrace Cool Shieldz**



- Filtration: 99+% filtration at .1 micron particles
- Anti fog eye shield (both sides)
- Fluid resistant face mask
- Attached Fluid Shield
- Fiberglass Free
- Latex Safe

Pleated, comfortable and fluid resistant protective face mask that breathes extremely well.

**Model No.**  
**11000-400 Ties (Surgical)**  
**11000-410 Ear Loop (Procedure)**

**One size**

**Embrace Protective Face Mask**



- Anti-Fog (Both Sides)
- Anti-Stat
- Optically clear
- Latex Safe







Soft elastic or adjustable Velcro strap  
**Embrace Splash Mask Face Shield**

**Model No.**  
**11000-000 Velcro Strap**  
**11000-001 Elastic Strap**





**One size**







# General & Utility Surgery Drapes

Drape	Description	Size	Adhesive	Use
	<b>D1000</b> Towel Drape	45 cm x 30 cm	45 cm x 3.8 cm	Prepping, tourniquet and instrument cover, squaring off, isolating sterile from non-sterile areas, fluid barrier, non-linting surface
	<b>D1000RL</b> <u>PREP-Drape®</u> Towel Drape on a roll	45 cm x 30 cm	45 cm x 3.8 cm	Prepping, tourniquet and instrument cover, isolating sterile from non-sterile areas, fluid barrier
	<b>D1003</b> Isolation Bag	50 cm x 50 cm		Isolate foot, hand, bowel, organs
	<b>D1010RL</b> <u>PREP-Drape®</u> Towel Drape on a roll	45 cm x 60 cm	45 cm x 3.8 cm	Prepping, tourniquet and instrument cover, isolating sterile from non-sterile areas, fluid barrier






Drape	Description	Size	Adhesive	Use
	<p><b>D1016</b> Fluid Drape</p>	50 cm x 48 cm	50 cm x 10 cm	CV, orthopedic, Neurologic, urologic and general
	<p><b>D1018</b> Instrument Pouch</p>	17 cm x 30 cm	(2) 17 cm x 5 cm	Neurologic, orthopedic, CV
	<p><b>D1040</b> Incise Drape</p>	60 cm x 35 cm	35 cm x 35 cm	Neurologic, orthopedic, back, CV, general
	<p><b>D1050</b> Incise Drape</p>	91 cm x 44 cm	59 cm x 44 cm	Neurologic, orthopedic, back, general, CV



Drape	Description	Size	Adhesive	Use
	<p><b>D1051</b> Incise Drape</p>	<p>91 cm x 83 cm</p>	<p>59 cm x 83 cm</p>	<p>CV, chest and back</p>
	<p><b>D1089</b> Utility Drape (2 Drapes per Pkg)</p>	<p>60 cm x 66 cm</p>		
	<p><b>D1092</b> Minor Procedure Drape</p>	<p>60 cm x 66 cm</p>	<p>6 cm circle</p>	<p>11 cm x 12 cm</p>
	<p><b>D2007</b> Clear Cover</p>	<p>91 cm x 111 cm</p>		



Drape	Description	Size	Aperture	Adhesive	Use
	<b>D2008</b> Clear Cover	111 cm x 182 cm			Back table cover, barrier sheet, set-up cover
	<b>D2147</b> Aperture Drape	137 cm x 139 cm	10 cm circle		ENT, general and utility
	<b>D2158</b> Instrument Cover	60 cm x 45 cm		60 cm x 1.2 cm	Cover instruments, mayo stand, monitor



# *Surgical Skin Marker*

## **TE-2751**

**Surgical Skin Marker, with regular tip**



**STERILE**



50/Box



## **6310**

**Surgical Skin Marker, TLS™ with regular tip, with Ruler**

**STERILE**

100/Box

