



GENERAL SURGICAL

TABLE OF CONTENTS

- ◆ Infection Control Glasses, Goggles, Face Masks, and Face Shields
- ◆ General & Utility Surgery Drapes
- ◆ Surgical Skin Marker

Infection Control Glasses, Goggles, Face Masks, & Face Shields



- Worn Alone
- Vented Sides
- Latex Safe
- Fits over prescription glasses
- Anti-Fog
- 99.9% UV protection



One size

Clear, lightweight plastic glasses provide wrap-around fluid protection and an extremely comfortable fit. Embrace Protective Glasses

Model No. 11000-300

Frame Clear



- UV400 Polycarbonate Lens
- Scratch Resistant Coating
- Lightweight Frames

Model No. 60S43CHS

Frame Clear



Complete peripheral protection with enhanced comfort.

SILO Wear Alone

One size



- Worn Alone
- Vented Sides
- Latex Safe
- Easily adjustable elastic strap
- Optical quality lenses for clarity of vision
- 99.9% UV protection

Sturdy, comfortable, snug-fitting for excellent fluid protection. Vented for additional comfort and reduced fogging.

Model No. 11000-328

Frame Clear

One size


Embrace Protective Goggles



5860-C Hollister Avenue
Santa Barbara, CA 93117 U.S.A

Tel: (805) 967 1153
Fax: (805) 967 1295

Home: www.tecfe.com
Email: tecfe@tecfe.com



- Worn Alone
- Wrap-around frame
- Secure & comfortable fit
- Non perforated lenses
- Anti-fog on both sides
- Top-down fluid protection ridge
- Fits over prescription glasses
- Latex Safe

The shield can be replaced in a matter of seconds providing extremely economical eye fluid protection.

Model No.
11000-500

Frame
Assorted colors

One size

Embrace Cool Shieldz



- Filtration: 99+% filtration at .1 micron particles
- Anti fog eye shield (both sides)
- Fluid resistant face mask
- Attached Fluid Shield
- Fiberglass Free
- Latex Safe

Pleated, comfortable and fluid resistant protective face mask that breathes extremely well.

Model No.
11000-400 Ties (Surgical)
11000-410 Ear Loop (Procedure)

One size

Embrace Protective Face Mask



- Anti-Fog (Both Sides)
- Anti-Stat
- Optically clear
- Latex Safe



Soft elastic or adjustable Velcro strap





Model No.
11000-000 Velcro Strap
11000-001 Elastic Strap

One size





Embrace Splash Mask Face Shield







General & Utility Surgery Drapes

| Drape | Description | Size | Adhesive | Use |
|---|---|---------------|----------------|---|
|  | D1000 Towel Drape | 45 cm x 30 cm | 45 cm x 3.8 cm | Prepping, tourniquet and instrument cover, squaring off, isolating sterile from non-sterile areas, fluid barrier, non-linting surface |
|  | D1000RL <u>PREP-Drape®</u> Towel Drape on a roll | 45 cm x 30 cm | 45 cm x 3.8 cm | Prepping, tourniquet and instrument cover, isolating sterile from non-sterile areas, fluid barrier |
|  | D1003 Isolation Bag | 50 cm x 50 cm | | Isolate foot, hand, bowel, organs |
|  | D1010RL <u>PREP-Drape®</u> Towel Drape on a roll | 45 cm x 60 cm | 45 cm x 3.8 cm | Prepping, tourniquet and instrument cover, isolating sterile from non-sterile areas, fluid barrier |






| Drape | Description | Size | Adhesive | Use |
|---|--|---------------|------------------|--|
|  | <p>D1016 Fluid Drape</p> | 50 cm x 48 cm | 50 cm x 10 cm | CV, orthopedic, Neurologic, urologic and general |
|  | <p>D1018 Instrument Pouch</p> | 17 cm x 30 cm | (2) 17 cm x 5 cm | Neurologic, orthopedic, CV |
|  | <p>D1040 Incise Drape</p> | 60 cm x 35 cm | 35 cm x 35 cm | Neurologic, orthopedic, back, CV, general |
|  | <p>D1050 Incise Drape</p> | 91 cm x 44 cm | 59 cm x 44 cm | Neurologic, orthopedic, back, general, CV |



| Drape | Description | Size | Adhesive | Use |
|---|--|-----------------------|----------------------|---------------------------|
|  | <p>D1051 Incise Drape</p> | <p>91 cm x 83 cm</p> | <p>59 cm x 83 cm</p> | <p>CV, chest and back</p> |
|  | <p>D1089 Utility Drape (2 Drapes per Pkg)</p> | <p>60 cm x 66 cm</p> | | |
|  | <p>D1092 Minor Procedure Drape</p> | <p>60 cm x 66 cm</p> | <p>6 cm circle</p> | <p>11 cm x 12 cm</p> |
|  | <p>D2007 Clear Cover</p> | <p>91 cm x 111 cm</p> | | |



| Drape | Description | Size | Aperture | Adhesive | Use |
|---|----------------------------------|-----------------|--------------|----------------|---|
|  | D2008 Clear Cover | 111 cm x 182 cm | | | Back table cover, barrier sheet, set-up cover |
|  | D2147 Aperture Drape | 137 cm x 139 cm | 10 cm circle | | ENT, general and utility |
|  | D2158 Instrument Cover | 60 cm x 45 cm | | 60 cm x 1.2 cm | Cover instruments, mayo stand, monitor |



Surgical Skin Marker & Corneal Marking Pad




TE-2751 Surgical Skin Marker

- ▶ Regular Tip
- ▶ Gentian Violet Ink
- ▶ 50 per box

STERILE



40490 Corneal Marking Pad



STERILE

- ▶ Gentian Violet Ink
- ▶ 10 per box



CE