

Nasal Sponges

AbsorbENT™ Nasal and Sinus Packing













We offer the most popular designs for nasal and sinus packing including Doyle Patient-Friendly Anatomic Nasal Packs. AbsorbENT™ Nasal and Sinus Packs consist of three-dimensional, interconnected open cells that are extremely soft and pliable. This interconnected structure produces packing that is exceptionally strong, highly absorbent and virtually lint and fiber free.

Soft, Absorbent, Lint and Fiber Free! AbsorbENT™ Nasal and Sinus Packs Feature

- Rounded edges for patient comfort during removal
- Even expansion for desired gentle mucosal pressure
- Rapid wicking of blood and fluid
- Individually sterile packaged

Nasal and Sinus Packing

AbsorbENT™ Nasal and Sinus Packs are compressed for easy placement. Upon exposure to fluids, the pack expands to six times its compressed size to fill the cavity and provide gentle pressure and support of tissues.

Product No.	Description	Packs/Box	Product Image
RH-2310-10	Epistaxis Pack with String, Large 10.0cm long x 1.5cm wide x 2.5cm thick	10	
RH-2315-20	Epistaxis Pack with String, Small 5.5cm long x 1.5cm wide x 2.5cm thick	20	
RH-7401-10	Nasal Packing 4.5cm long x 1.5cm wide x 2.0 cm high	10	
RH-7404-10	Nasal Packing 8.0cm long x 1.5cm wide x 2.0cm high	10	
RH-7405-10	Nasal Packing with Airway and String 8.0cm long x 1.5cm wide x 2.0cm thick	10	
RH-7408-10	Nasal Packing with Airway and String 4.5cm long x 1.5cm wide x 2.0cm thick	10	
RH-7409-10	Nasal Packing, Anatomical with String 8.0cm long x 1.5cm wide x 3.0cm thick	10	
RH-7410-10	Nasal Packing, Anatomical with Airway and String 8.0cm long x 1.5cm wide x 3.0cm thick	10	
RH-7406-10	Nasal Packing, Thin 8.0cm long x 1.0cm wide x 2.0cm thick	10	
RH-7411-10	Nasal Packing, Thin Anatomical 8.0cm long x 1.0cm wide x 3.0cm thick	10	
RH-7530-10	Sinus Packing with String 3.5cm long x 0.6cm wide x 1.2cm thick	10*	
RH-7532-10	Sinus Packing with String 3.5cm long x 0.9cm wide x 1.2cm thick	10*	

Doyle Patient Friendly™ Nasal Packs



Tecfen Corporation

5860-C Hollister Avenue
Santa Barbara, CA 93117 U.S.A



Tel: (805) 967 1153
Fax: (805) 967 1295



Home: www.tecfe.com
Email: tecfe@tecfe.com

We offer the most popular designs for nasal and sinus packing including Doyle Patient-Friendly Anatomic Nasal Packs. AbsorbENT™ Nasal and Sinus Packs consist of three-dimensional, interconnected open cells that are extremely soft and pliable. This interconnected structure produces packing that is exceptionally strong, highly absorbent and virtually lint and fiber free.

Soft, Absorbent, Lint and Fiber Free! AbsorbENT™ Nasal and Sinus Packs Feature



- Rounded edges for quick and easy insertion
- Even expansion for desired gentle mucosal pressure
- Rapid wicking of blood and fluid
- Individually sterile packaged

Doyle Patient Friendly™ Nasal Packs*

Doyle Patient-Friendly Nasal Packs are sterile, anatomically shaped, compressed AbsorbENT™ Sponges. The pack is precut from back to front except for a web at the distal end. The pack is inserted and functions like a standard anatomical nasal pack. On removal, the pack is pulled slightly out of the nasal cavity to expose the uncut web. The web is then cut through with

scissors, which separates the pack into two parts. The two parts are removed one at a time, greatly reducing patient discomfort.



Inserted into nasal fossa



Pull out slightly and cut through remaining web



Remove two parts one at a time

*Developed in cooperation with Donald Doyle, MD, FACS

Product No.	Description	Packs/Box	Product Image
RH-7412-10	Doyle Patient Friendly Anatomical Nasal Packing with String 8.0cm long x 1.5cm wide x 3.0cm thick	10	
RH-7413-10	Doyle Patient Friendly Anatomical Nasal Packing with Airway Tube 8.0cm long x 1.5cm wide x 3.0cm thick	10	
RH-7414-10	Doyle Patient Friendly Anatomical Nasal Packing with String, Thin 8.0cm long x 1.0cm wide x 3.0cm thick	10	



Tecfen Corporation

5860-C Hollister Avenue
Santa Barbara, CA 93117 U.S.A



Tel: (805) 967 1153
Fax: (805) 967 1295



Home: www.tecfe.com
Email: tecfe@tecfe.com