

Dual Eyeshields, Goggles, Single Eye Shields & Moisture Chamber



5818D CLEAR



5819D SMOKE
with UV Protection

Encore® Dual Eyeshield

Designed to be worn at night by post LASIK and cataract patients for maximum comfort & protection

- ▶ Increases all-around optical visibility
- ▶ Adjustable to fit all head sizes
- ▶ Latex free
- ▶ Made of polycarbonate and durable single-cell foam



5820GS SMOKE
with UV Protection



5821GC CLEAR
with UV Protection

Encore® Goggle Eyeshield

Polycarbonate lenses provide maximum protection and comfort for dry eye & post operative patients

- ▶ Adjustable head strap
- ▶ Latex free
- ▶ Hinged nose bridge provides added comfort
- ▶ Durable foam on shield is slip proof & comfortable



5702CS CLEAR



5703CS SMOKE
with UV Protection

Encore® Single Eyeshield

Polycarbonate allows for maximum comfort and protection for post-operative patients

- ▶ Center design provides maximum optical clarity
- ▶ Specially designed ventilation prevents fogging



TE-501040 CLEAR,
STERILE

- ▶ Individually packaged
- ▶ 50 pouches/bag

9401 CLEAR,
Non-STERILE

- ▶ 100/bag



5704C



5824MS

Encore® 1 and 2 Piece Eyeshield

Polycarbonate lenses provide maximum protection and comfort for post-operative patients

- ▶ Adjustable head strap
- ▶ Durable foam is slip proof and comfortable
- ▶ Latex free
- ▶ Ventilated to reduce fogging



5822MS CLEAR

Encore® Moisture Chamber

Designed to be worn at night to retain moisture for dry eye

- ▶ Adjustable to fit all head sizes
- ▶ Made of polycarbonate and durable single-cell foam
- ▶ Latex free

