

Eyeshields



5818D CLEAR



5819D SMOKE *with UV Protection*

Encore® Dual Eyeshield

- ▶ Designed to be worn at night by post LASIK and cataract patients for maximum comfort & protection
- ▶ Increases all-around optical visibility
- ▶ Adjustable to fit all head sizes
- ▶ Latex free
- ▶ Made of polycarbonate and durable single-cell foam



5821GC CLEAR
with UV Protection



5820GS SMOKE
TINT *with UV Protection*

Encore® Goggle Eyeshield

- ▶ Polycarbonate lenses provide maximum protection and comfort for dry eye & post operative patients
- ▶ Adjustable head strap
- ▶ Latex free
- ▶ Hinged nose bridge provides added comfort
- ▶ Durable foam on shield is slip proof & comfortable

Encore® Single Eyeshield

- ▶ Polycarbonate allows for maximum comfort and protection for post-operative patients
- ▶ Center design provides maximum optical clarity
- ▶ Specially designed ventilation prevents fogging



5702CS
CLEAR




5703CS
SMOKE TINT
with UV Protection



RTM Eyeshields

- 529** Newborn (18.4 cm)
- 530** Premature (13.3 cm)
- 531** Small Premature (8.2 cm)



- ▶ Soft knit fabric laminated over dark gray filter to assure complete opacity
- ▶ Covers only eyes eliminating interference with cephalic IV's and any risk of airway occlusion
- ▶ Patented strapless phototherapy protection. No straps to adjust
- ▶ Tabs have sensitivity test 3M™ hypoallergenic adhesive that is safe on infant's skin

