













# General & Utility Surgery Drapes

Drape	Description	Size	Adhesive	Use
	<b>D1000</b> Towel Drape	45 cm x 30 cm	45 cm x 3.8 cm	Prepping, tourniquet and instrument cover, squaring off, isolating sterile from non-sterile areas, fluid barrier, non-linting surface
	<b>D1000RL</b> <u>PREP-Drape®</u> Towel Drape on a roll	45 cm x 30 cm	45 cm x 3.8 cm	Prepping, tourniquet and instrument cover, isolating sterile from non-sterile areas, fluid barrier
	<b>D1003</b> Isolation Bag	50 cm x 50 cm		Isolate foot, hand, bowel, organs
	<b>D1010RL</b> <u>PREP-Drape®</u> Towel Drape on a roll	45 cm x 60 cm	45 cm x 3.8 cm	Prepping, tourniquet and instrument cover, isolating sterile from non-sterile areas, fluid barrier






Drape	Description	Size	Adhesive	Use
	<p><b>D1016</b> Fluid Drape</p>	50 cm x 48 cm	50 cm x 10 cm	CV, orthopedic, Neurologic, urologic and general
	<p><b>D1018</b> Instrument Pouch</p>	17 cm x 30 cm	(2) 17 cm x 5 cm	Neurologic, orthopedic, CV
	<p><b>D1040</b> Incise Drape</p>	60 cm x 35 cm	35 cm x 35 cm	Neurologic, orthopedic, back, CV, general
	<p><b>D1050</b> Incise Drape</p>	91 cm x 44 cm	59 cm x 44 cm	Neurologic, orthopedic, back, general, CV



Drape	Description	Size	Adhesive	Use
	<p><b>D1051</b> Incise Drape</p>	<p>91 cm x 83 cm</p>	<p>59 cm x 83 cm</p>	<p>CV, chest and back</p>
	<p><b>D1089</b> Utility Drape (2 Drapes per Pkg)</p>	<p>60 cm x 66 cm</p>		
	<p><b>D1092</b> Minor Procedure Drape</p>	<p>60 cm x 66 cm</p>	<p>6 cm circle</p>	<p>11 cm x 12 cm</p>
	<p><b>D2007</b> Clear Cover</p>	<p>91 cm x 111 cm</p>		



Drape	Description	Size	Aperture	Adhesive	Use
	<b>D2008</b> Clear Cover	111 cm x 182 cm			Back table cover, barrier sheet, set-up cover
	<b>D2147</b> Aperture Drape	137 cm x 139 cm	10 cm circle		ENT, general and utility
	<b>D2158</b> Instrument Cover	60 cm x 45 cm		60 cm x 1.2 cm	Cover instruments, mayo stand, monitor

