

PVA & Cellulose Sponge Eye Products

Sterile **PVA** Eye Spears

Lint-free, highly absorbent, specifically designed for LASIK & other delicate surgical procedures

(Also Available in Non-STERILE)

Item No.	Description
TE-091B	5 spears/foil ♦ 30 foils/box = 150 spears/box
TE-091D	5 spears/foil ♦ 100 foils/bag = 500 spears/bag

Sterile **Cellulose** Eye Spears

Lint-free, highly absorbent, specifically designed for LASIK & other delicate surgical procedures

(Also Available in Non-STERILE)

40415 10 spears/glassine ♦ 1 glassine/pack ♦ 20 packs/box = 200 spears/box



Lasik Disk

40804 UltraCell Sterile Lasik Disk w/Dark Blue Handle (9mm diameter 2mm thickness), 2/Pack, 20 Packs/Box



Lasik Shield

40820 Lasik Shield Sterile (7mm flat spot), 1/Pack, 20/Box



40821

Lasik Shield Sterile (6mm flat spot), 1/Pack, 20/Box

40822

Lasik Shield Sterile (4 and 8mm flat spot), 1/Pack, 20/Box



Lasik Drain

40830

Lasik Drain Sterile (11.5mm opening) 42mm x 17mm, 1/Pack, 20/Box

40831

Lasik Drain Sterile (10.5mm opening) 40mm x 14mm, 1/Pack, 20/Box



Corneal Light Shield

40420

Corneal Light Shield Sterile (7mm), 1/Pack, 20/Box

40421

Corneal Light Shield Sterile (8mm), 1/Pack, 20/Box

40422

Corneal Light Shield Sterile (9mm), 1/Pack, 20/Box



Eye Drain & Wick

40430

Eye Drain Sterile (80cc capacity), 1/Pack, 20/Box

40431

Eye Wick Sterile (4mm x 170mm), 2/Pack, 10/Box

40435

Eye Drain Sterile (400cc capacity), 1/Pack, 10/Box



Ophthalmic Ring

40440

Bloomberg Ophthalmic Ring Sterile (1.85cm), 1/Pack, 20/Box





Eye Kit

40450

Eye Kit Sterile, 1 wipe/wick & 1 Light Shield (7mm),
1/Pack, 20/Box



Cataract Kit

40451

Cataract Kit Sterile, (5 Cell. Spears, 1 Corn. Light Shield
(7mm), 1 Inst. Wipe (2mm), 1 Eye Drain (80cc), 1/Pack,
20/Box



Eye Pledgets

40825

Eye Pledgets Sterile (2mm x 6mm x 1mm thick),
10/Pack, 20/Box

Instrument Wipe



90910

Instrument Wipe Sterile (7.6cm x 7.6 cm x 1mm),
1/Pack, 20/Box

