













# Ophthalmic Surgery Drapes

Drape	Description	Size	Aperture	Adhesive	Procedure
	<b>D1020</b> Aperture Drape	40 cm x 40 cm	6 cm circle	12 cm x 11 cm	Cataract, corneal, refractive, lid/plastic, minor
	<b>D1021</b> Incise Drape	40 cm x 40 cm		12 cm x 11 cm	Cataract, refractive, vitreoretinal, corneal
	<b>D1022</b> Refractive Drape	40 cm x 27 cm	3.8 cm x 2.5 cm	7.6 cm x 7.6 cm	Lasik, refractive, corneal, lid/plastic
	<b>D1024</b> Ophthalmic Aperture Pouch	48 cm x 38 cm	5 cm x 3 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal







Drape	Description	Size	Aperture	Adhesive	Procedure
	<b>D1025</b> Ophthalmic Incise Pouch	48 cm x 38 cm		12 cm x 11 cm (2) 2.5 cm x 5 cm TABS	Cataract, vitreoretinal, corneal
	<b>D1030</b> Aperture Drape	83 cm x 129 cm	7 cm x 5 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal, lid/plastic,
	<b>D1031</b> Ophthalmic Aperture Drape Pouch	121 cm x 129 cm	7 cm x 5 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal
	<b>D1035</b> Incise Drape	16 cm x 19 cm		10 cm x 19 cm	Cataract, vitreoretinal, corneal, lid/plastic, minor



Drape	Description	Size	Aperture	Adhesive	Procedure	
	<b>D1060</b> Ophthalmic Incise Drape	121 cm x 129 cm			12 cm x 11 cm	Cataract, vitreoretinal, corneal
	<b>D1061</b> Ophthalmic Incise Drape <b>Pouch</b>	121 cm x 129 cm			12 cm x 11 cm	Cataract, vitreoretinal
	<b>D1062</b> Ophthalmic Incise Drape <b>Pouch</b>	121 cm x 139 cm			12 cm x 11 cm	Cataract, vitreoretinal
	<b>D1064</b> Ophthalmic Aperture Sheet <b>Pouch</b>	139 cm x 177 cm	5 cm x 3.8 cm		11 cm x 12 cm	Cataract, vitreoretinal,





Drape	Description	Size	Aperture	Adhesive	Procedure	
	<b>D1065</b> Ophthalmic Incise Drape <b>Pouch</b>	139 cm x 177 cm			11 cm x 12 cm	Cataract, vitreoretinal,
	<b>D1122</b> Refractive Drape	40 cm x 27 cm	3.8 cm x 2.54 cm		7.6 cm x 7.6 cm	Lasik, refractive, corneal, lid/plastic
	<b>D1164</b> Ophthalmic Aperture Sheet <b>Pouch</b>	139 cm x 292 cm	5 cm x 3.8 cm		11 cm x 12 cm	Cataract, vitreoretinal,
	<b>D1165</b> Ophthalmic Incise Drape <b>Pouch</b>	139 cm x 292 cm			11 cm x 12 cm	Cataract, vitreoretinal,



Drape	Description	Size	Aperture	Adhesive	Procedure
	<b>D2019</b> Ophthalmic Incise Sheet	106 cm x 121 cm		7 cm circle	Cataract, vitreoretinal, corneal
	<b>D2020</b> Ophthalmic Aperture Drape	106 cm x 121 cm	7 cm x 5 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal
	<b>D2021</b> Irrigation Pouch	21 cm x 27 cm		21 cm x 3.8 cm	Cataract, vitreoretinal, corneal, ER, minor
	<b>D2043</b> Strabismus Drape	101 cm x 121 cm	13 cm x 3.8 cm	17 cm x 11 cm	Strabismus, lid/plastic, corneal



Drape	Description	Size	Aperture	Adhesive	Procedure
	<p><b>D2114</b> Irrigation Pouch</p>	<p>30 cm x 27 cm</p>		<p>30 cm x 3.8 cm</p>	<p>Cataract, vitreoretinal, corneal, ER, minor</p>
	<p><b>D2258</b> Ophthalmic Incise Sheet Pouch</p>	<p>129 cm x 190 cm</p>		<p>12 cm x 11 cm</p>	<p>Cataract, vitreoretinal, corneal</p>

