








Ophthalmic Surgery Drapes

Drape	Description	Size	Aperture	Adhesive	Procedure
	D1020 Aperture Drape	40 cm x 40 cm	6 cm circle	12 cm x 11 cm	Cataract, corneal, refractive, lid/plastic, minor
	D1021 Incise Drape	40 cm x 40 cm	12 cm x 11 cm		Cataract, refractive, vitreoretinal, corneal
	D1022 Refractive Drape	40 cm x 27 cm	3.8 cm x 2.5 cm	7.6 cm x 7.6 cm	Lasik, refractive, corneal, lid/plastic
	D1024 Ophthalmic Aperture Pouch	48 cm x 38 cm	5 cm x 3 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal



Drape	Description	Size	Aperture	Adhesive	Procedure
	D1025 Ophthalmic Incise Pouch	48 cm x 38 cm		12 cm x 11 cm (2) 2.5 cm x 5 cm TABS	Cataract, vitreoretinal, corneal
	D1030 Aperture Drape	83 cm x 129 cm	7 cm x 5 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal, lid/plastic,
	D1031 Ophthalmic Aperture Drape POUCH: 30 cm x 27 cm	121 cm x 129 cm	7 cm x 5 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal
	D1035 Incise Drape	16 cm x 19 cm		10 cm x 19 cm	Cataract, vitreoretinal, corneal, lid/plastic, minor



Drape	Description	Size	Aperture	Adhesive	Procedure	
	D1060 Ophthalmic Incise Drape	121 cm x 129 cm			12 cm x 11 cm	Cataract, vitreoretinal, corneal
	D1061 Ophthalmic Incise Drape POUCH: 30 cm x 27 cm	121 cm x 129 cm			12 cm x 11 cm	Cataract, vitreoretinal
	D1062 Ophthalmic Incise Drape POUCH: 30 cm x 27 cm	121 cm x 139 cm			12 cm x 11 cm	Cataract, vitreoretinal
	D1064 Ophthalmic Aperture Sheet POUCH: 55 cm x 29 cm	139 cm x 177 cm	5 cm x 3.8 cm		11 cm x 12 cm	Cataract, vitreoretinal,






Drape	Description	Size	Aperture	Adhesive	Procedure
	D1065 Ophthalmic Incise Drape POUCH: 55 cm x 29 cm	139 cm x 177 cm		11 cm x 12 cm	Cataract, vitreoretinal,
	D1122 Refractive Drape	40 cm x 27 cm	3.8 cm x 2.54 cm	7.6 cm x 7.6 cm	Lasik, refractive, corneal, lid/plastic
	D1164 Ophthalmic Aperture Sheet POUCH: 55 cm x 29 cm	139 cm x 292 cm	5 cm x 3.8 cm	11 cm x 12 cm	Cataract, vitreoretinal,
	D1165 Ophthalmic Incise Drape POUCH: 55 cm x 29 cm	139 cm x 292 cm		11 cm x 12 cm	Cataract, vitreoretinal,



Drape	Description	Size	Aperture	Adhesive	Procedure
	D2019 Ophthalmic Incise Sheet	106 cm x 121 cm		7 cm circle	Cataract, vitreoretinal, corneal
	D2020 Ophthalmic Aperture Drape	106 cm x 121 cm	7 cm x 5 cm	12 cm x 11 cm	Cataract, vitreoretinal, corneal
	D2021 Irrigation Pouch	21 cm x 27 cm		21 cm x 3.8 cm	Cataract, vitreoretinal, corneal, ER, minor
	D2043 Strabismus Drape	101 cm x 121 cm	13 cm x 3.8 cm	17 cm x 11 cm	Strabismus, lid/plastic, corneal



Drape	Description	Size	Aperture	Adhesive	Procedure
	D2114 Irrigation Pouch	30 cm x 27 cm		30 cm x 3.8 cm	Cataract, vitreoretinal, corneal, ER, minor
	D2258 Ophthalmic Incise Sheet POUCH: 55 cm x 29 cm	129 cm x 190 cm		12 cm x 11 cm	Cataract, vitreoretinal, corneal



PD310
Plastic Pouch Drape
60cm x 40cm
with incise



PD440 Half Body
Non-Woven Drape
140cm x 127cm
with incise & pouch



PD420 Half Body
Non-Woven Drape
140cm x 127cm with 6cm
oval aperture & pouch



PD620 Full Body
Non-Woven Drape
165cm x 280cm
with incise & 2
pouches

