

# ODYSSEY Punctal Plugs & Intracanalicular Plugs

Patented products for the treatment of dry eye syndrome  
for Permanent, Temporary, and Partial Occlusion

## PERMANENT OCCLUSION

### Parasol® Sterile Pre-Loaded Punctal Plugs (1 pair/box)

A14-200	X-Small (0.25mm-0.35mm)
A14-201	Small (0.35mm-0.65mm)
A14-202	Medium (0.6mm-0.85mm)
A14-203	Large (0.9mm+)

### Parasol® Sterile QuickLoad™ Punctal Plugs (8 pair/box)


A14-201/QL	Small (0.35mm-0.65mm)
A14-202/QL	Medium (0.6mm-0.85mm)
A14-203/QL	Large (0.9mm+)
A14-205/QL	Combo Pack (4 Pr Small, 4 Pr Medium)

### Parasol® Non-Sterile Bulk Punctal Plugs (10 pair/box)

A14-201BN	Small (0.35mm-0.65mm)
A14-202BN	Medium (0.6mm-0.85mm)
A14-203BN	Large (0.9mm+)
A14-209BN	Combo Pack (5 Pr Small, 5 Pr Medium)

### Parasol® Non-Sterile Bulk Punctal Plugs (5 pair/box)

A20-203BN	Large (0.9mm+)
-----------	----------------



## PARTIAL OCCLUSION

The Micro Flow™ punctal occluder is designed to provide a simple, reversible means of partial occlusion by limiting the drainage of tears into the lacrimal duct.


### Micro Flow™ Sterile Pre-Loaded Punctal Plugs (1 pair/box)

A14-401	Small (0.4mm-0.55mm)
A14-402	Medium (0.55mm-0.7mm)
A14-403	Large (0.7mm-0.85mm)



### Micro Flow™ Non-Sterile Punctal Plugs (5 pair/box)

A20-401/BN	Small (0.4mm-0.55mm)
A20-402/BN	Medium (0.55mm-0.7mm)
A20-403/BN	Large (0.7mm-0.85mm)



## TEMPORARY OCCLUSION

### Temporary Collagen Implants (60 implants/box)

Provide occlusion for 7-10 days for diagnostic use

A12-102	0.2mm X 2.0mm
A12-103	0.3mm X 2.0mm
A12-104	0.4mm X 2.0mm
A12-200	Variety Pack : (18 implants 0.2mm, 24 implants 0.3mm, 18 implants 0.4mm)

### EXTEND™ Synthetic Absorbable Implants (10 pair/box)

Provide occlusion for up to 3 months

A13-202	0.2mm x 2.0mm
A13-203	0.3mm x 2.0mm
A13-204	0.4mm x 2.0mm
A13-205	0.5mm x 2.0mm
A13-200	Variety Pack: (3 pair 0.3mm, 5 pair 0.4mm, 2 pair 0.5mm)



#### Insertion Instrument

A14-300 1 pair per box



#### Forceps

A14-301 1 pair Per box



5860-C Hollister Avenue  
Santa Barbara, CA 93117 U.S.A.

Telephone: 805 967-1153  
FAX: 805 967-1295

Website: [www.tecfe.com](http://www.tecfe.com)  
email: [tecfe@tecfe.com](mailto:tecfe@tecfe.com)