

# *Tecfen Disposable Blanket*

## TE-105 TECFEN DISPOSABLE BLANKET



The Disposable Blanket is highly absorbent, soft and pliable blanket made of 2-ply bleached non-wet strength paper tissue laminated onto a thick yellow low density polyethylene film. The blanket provides excellent protection, padding, and cushioning. In emergency situations, it is used because of its high absorption of body fluids.

### Uses:

- Non-sterile emergency medical use
- Protective underlay for many different uses
- Emergency Blanket for body heat retention
- Assists in prevention of post-accident shock
- Emergency Tarp for weather protection
- Reversible for liquid absorption or water repellence
- Soft pliable material
- Light weight & compact packaging

### Specifications:

Color: Yellow  
Size: 137cm x 203cm  
Weight: 180 grams



### **Tecfen Corporation**

5860-C Hollister Avenue  
Santa Barbara, CA 93117 U.S.A

Tel: (805) 967 1153  
Fax: (805) 967 1295



**Home:** [www.tecfe.com](http://www.tecfe.com)  
Email: [tecfe@tecfe.com](mailto:tecfe@tecfe.com)