

# Infection Control Glasses, Goggles, Face Masks, & Face Shields



- Worn Alone
- Vented Sides
- Latex Safe
- Fits over prescription glasses
- Anti-Fog
- 99.9% UV protection



One size

Clear, lightweight plastic glasses provide wrap-around fluid protection and an extremely comfortable fit. Embrace Protective Glasses

Model No. 11000-300

Frame Clear



- UV400 Polycarbonate Lens
- Scratch Resistant Coating
- Lightweight Frames

Model No. 60S43CHS

Frame Clear



Complete peripheral protection with enhanced comfort.

SILO Wear Alone

One size



- Worn Alone
- Vented Sides
- Latex Safe
- Easily adjustable elastic strap
- Optical quality lenses for clarity of vision
- 99.9% UV protection

Sturdy, comfortable, snug-fitting for excellent fluid protection. Vented for additional comfort and reduced fogging.

Model No. 11000-328

Frame Clear

One size


Embrace Protective Goggles



5860-C Hollister Avenue  
Santa Barbara, CA 93117 U.S.A

Tel: (805) 967 1153  
Fax: (805) 967 1295

Home: [www.tecfe.com](http://www.tecfe.com)  
Email: [tecfe@tecfe.com](mailto:tecfe@tecfe.com)



- Worn Alone
- Wrap-around frame
- Secure & comfortable fit
- Non perforated lenses
- Anti-fog on both sides
- Top-down fluid protection ridge
- Fits over prescription glasses
- Latex Safe

The shield can be replaced in a matter of seconds providing extremely economical eye fluid protection.

**Model No.**  
**11000-500**

**Frame**  
**Assorted colors**

**One size**

### Embrace Cool Shieldz



- Filtration: 99+% filtration at .1 micron particles
- Anti fog eye shield (both sides)
- Fluid resistant face mask
- Attached Fluid Shield
- Fiberglass Free
- Latex Safe

Pleated, comfortable and fluid resistant protective face mask that breathes extremely well.

**Model No.**  
**11000-400 Ties (Surgical)**  
**11000-410 Ear Loop (Procedure)**

**One size**

### Embrace Protective Face Mask



- Anti-Fog (Both Sides)
- Anti-Stat
- Optically clear
- Latex Safe



Soft elastic or adjustable Velcro strap

**Model No.**  
**11000-000 Velcro Strap**  
**11000-001 Elastic Strap**

**One size**

### Embrace Splash Mask Face Shield

